



Helps you
recommend
**the best course
of treatment**
for your patient.

Fast Confident and
Completely Personalized

TECHNICAL SPECIFICATIONS


Power Supply	Input: 100-240V, 50-60Hz	Output: 12V, 1-2A
Max. Dimensions	Handset: 15x6x7 (cm)	Dock: 19x21x10 (cm)
Weight	Handset: 150g approx.	Dock: 661g approx.
Measurement method	Specialized level index (0-99) for each measurement category	
Image Modes	Hydration, Pores, Spots/Pigmentation, Wrinkles, Acne, Keratin & Sebum	
Specialized Lens for Skin Measurement	5 Mega Pixel Image Sensor Various image modes with more than 20 LEDs Real-time image transfer and easy to handle	
Support Protocol	USB 2.0 • Wi-Fi IEEE 802.11 b/g/n • NFC	
Certification		

Solely Distributed by: _____

DUBAI Medica Trading LLC Tel: + 971 4 4370 669 Fax: + 971 4 4370 668 info@medicagroup.com	ABU DHABI Medica Trading LLC Tel: + 971 2 6581 619 Fax: + 971 2 6581 612 info@medicagroup.com	KSA Medica - RIC Tel: + 966 11 4654 113 Fax: + 966 11 4630 135 info@medicagroup.com	LEBANON Medica Pharmaceuticals s.a.l. Tel: + 961 4 710 973 Fax: + 961 4 710 977 info-me@medicagroup.com	SINGAPORE Medica Singapore Pte Ltd Tel: + 65 67 33 96 89 Fax: + 65 67 33 96 91 info-sg@medicagroup.com
------------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------

Skin & Hair
Analysis System

Fast & Completely Personalized Customer Service





A new professional system for measurement, cosmetic treatment and recommendation of skin and hair products.

SRS DermAnalysis is the leading and most advanced range of measuring equipment for skin and hair health condition, which includes fast and completely personalized customer service, focused on the customer's needs.

SRS DermAnalysis system tells you the course of treatment and products that may work best for each customer based on their analysis results.

It is a single device with built-in lenses and does not require any extra accessories which makes the analysis process more convenient and as fast as possible.

- ✓ User friendly
- ✓ Easy to use camera with polarized lens
- ✓ Powerful customer management software
- ✓ Marketing tool to help you promote and sell dermocosmetic products
- ✓ Several applications for skin and hair
- ✓ Remote operation through Android tablets with the SRS Skin - DermAnalysis Software

PRECISE & FAST

A triple-mode image processing lens is used to automatically conduct a complete analysis of the desired area.

Images can be compared to analyse specific details of the skin and hair's condition

The results are shown on-screen with graphics that customers can easily interpret.



- Real-time image transfer via WiFi.
- E-mail sent with personalized diagnosis results.



TAILORED TREATMENT

Shows the treatment needed and the products that may work best for each customer based on the analysis results and their age segment.

The results of each analysis are saved in the system to facilitate personal follow-up

Professionals can customize the products recommended by the system in line with their offers, promotions or sales agreements.

SRS DermAnalysis automatically carries out a complete diagnosis of the main parameters

FOR SKIN

HYDRATION

PORES

SPOTS

SKIN CLARITY

SKIN TONE

WRINKLES

ACNE

KERATIN

FOR HAIR

LOSS TYPE

DENSITY

SENSITIVITY

THICKNESS

DEAD SKIN CELLS



High Resolution Wireless Diagnostic System

CAMERA

Wireless video transmission in real time, with high magnification images.

image sensor	5MEGA Pixel CMOS with Optic Zoom		
Connection	Wifi, Wifi direct		
Weight	aprox 150g		
Battery	1050 mAh (Li-Polymer), 60-90 minutes of continuous use. Fast charging connection through the base (12V)		
Image Mode	Skin	VSL Mode	Pores - Wrinkles
		XPL Mode	Spots - Keratin
		UVL Mode	Acne
Integrated Hydration Sensor			
	Hair	VSL Mode	Density
		XPL Mode	Sensitivity
			Dead Skin Cells
			Thickness



Self Sterilizing for peace of mind!

